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Chastru Biotech Limited.
Add: 2/F, Yau Tak Building, 167 Lockhart Road, Wanchai, Hong Kong, China.

SUNGO Europe B.V.
Add: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands Contact: Yana Zhang (Ms.), Tel/Fax: +31 (0) 2021 11106

E-Mail: yan.zhang@sungoglobal.com

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Chargennummer

Verfallsdatum

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Hersteller

Symbol Index

Professionelle In-Vitro-Diagnostik

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Chastru COVID-19 Antigen Rapid Test 4 in 1 ENGLISH

For in vitro diagnostic use only

[INTENDED USE]

The Chastru COVID-19 Antigen Rapid Test 4in1 is a lateral flow immunoassay for a qualitative detection of SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swabs, oropharyngeal swabs, nasal swabs (anterior nasal) and saliva samples in individuals, who are under the suspicion on a SARS-CoV-2 infection.

The test results relate to the identification of the SARS-CoV-2 nucleocapsid antigen. The antigen is generally in detectable in samples taken during the acute phase of the infection. Positive results indicate the presence of viral antigens. A clinical correlation with the medical history and other diagnostic information is necessary, in order to determine the infection status. Positive results do not rule out a bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Negative results do not exclude a SARS-CoV-2-infection and should not be used as the sole basis for treatment or for patient management decisions, including decisions for infection control. Negative results should be considered in relation to a patient's recent exposures, history and presence of clinical signs and symptoms consistent with COVID-19 and, where appropriate, confirmed by an RT-PCR laboratory test. The COVID-19 Antigen Rapid Test Kit is intended for the use by trained clinical personnel, specially prepared and trained for in vitro diagnostic procedures.

[SUMMARY]

The new coronavirus (SARS-CoV-2) belongs to the genus β . COVID-19 is an acute respiratory infection disease caused by the virus. Currently, patients infected with the new Coronavirus are the main source of infection. Asymptotically infected people can also be an infectious source. According to the current epidemiological investigation, the incubation period is 1 to 14 days, usually 3 to 7 days. The symptoms are fever, tiredness and dry cough. In some cases, a blocked or runny nose, sore throat, muscle pain and diarrhea occur.

[PRINCIPLE]

The COVID-19 Antigen Rapid Test Kit is a chromatographic immunoassay kit based on the principle of the double antibody technology. The COVID-19 Antigen Rapid Test Kit is used to detect nucleocapsid antigens from SARS-CoV-2 in nasopharyngeal swabs, oropharyngeal swabs, nasal swabs (anterior nasal) and saliva samples from patients suspected of having COVID-19. During the test, the test membrane strip is moistened by the sample material and the capillary action. The SARS-CoV-2 antigens that may be present in the sample bind to the antibody conjugates. The immune complex is then caught by monoclonal SARS-CoV-2 nucleocapsid protein antibodies on the pre-coated membrane. As a consequence, a colored line within the test line range appears; this indicates a positive result. In the absence of SARS-CoV-2 antigens, no colored line is formed in the test line area, which indicates a negative result. For procedural control, a colored line always appears in the control line area. This indicates that the correct sample volume has been added and a sufficient membrane wetting by capillary action has taken place.

[PRECAUTIONS]

- Only for in vitro diagnostics by professionals, such as health professionals and nursing professionals.
- Do not use this product as the sole basis for diagnosing or excluding SARS-CoV-2 infection.
- Do not use after expiration date.
- Read all information in this package before you start the test.
- The test kit should remain in the sealed bag until use.
- Treat all samples as potentially infectious. Observe proven precautions against microbiological hazards throughout the test.
- Dispose all used individual components of the rapid test in accordance with local regulations.

[MATERIALS SUPPLIED]

The test kit contains a membrane strip, which is coated with a monoclonal anti-SARS-CoV-2 nucleocapsid protein antibody on the T Test line, as well as a coloring material pad, which contains colloidal gold and is coupled to a monoclonal SARS-CoV-2-nucleocapsid protein antibody. The number of tests contained in the test pack was printed on the test package packaging.

Materials supplied

- Test kits - individually sealed
- Sterilized sampling swabs
- Extraction tube for sample extraction, filled with extraction solutions
- Dropper caps for sample tubes
- Sample tube holder
- Instructions for use

Materials required but not provided: Timer

[CONSERVATION AND STABILITY]


- Store at room temperature or refrigerated between 4°C and 30°C in the package and sealed bag. The test can be kept until the expiry date stated on the sealed bag.
- After opening the bag, the test should be used within one hour. Prolonged contact with hot and humid environment leads to product deterioration.
- The lot number (batch number) and the expiration date were printed on the packaging of the sealed test kits.

[SAMPLE COLLECTION AND HANDLING]

Samples, which are taken in the early symptom phase, contain the highest virus titers; samples, which are taken after five days with symptoms, tend to give negative results compared to an RT-PCR test. An insufficient sampling, an inappropriate handling and/or inappropriate transport of the samples can lead to a false negative result. Professional training in sampling is therefore strongly recommended, since sample quality is important for the preparation of accurate test results.

Sample collection

Nasopharyngeal swab sample

Carefully insert the flexible sampling swab into the 70° tilted head parallel to the palate (not upward) through the nostril until resistance occurs or the distance of the inserted swab equals that from the ear to the patient's nostril, indicating contact with the nasopharyngeal zone. The swab should reach a depth, which corresponds to the distance between the nostrils and the outside opening of the ear. Rub and roll the swab carefully. Leave the swab in place for a few seconds to absorb secretions. Turn the swab while slowly removing it.


Samples can be taken from both nostrils with the same swab. However, it is not necessary to collect samples from both nostrils when the tip of the swab is saturated with liquid from the first sample collection. If a curved nasal septum or a blockage causes difficulties in taking the sample out of a nostril, use the same swab to take the sample out of the other nostril.

Oropharyngeal swab sample

Insert the sampling swab into the posterior pharynx and tonsil area. Rub the swab over both tonsils and posterior oropharynx and avoid touching the tongue, teeth and gums.

Nasal swab sample - anterior nasal

Carefully insert the flexible sampling swab approx. 2.5 cm into a nostril until a slight resistance is felt. Gently rub and roll the swab along the inner wall of the nose for at least 15 seconds, taking up as much sample material as possible. While slowly removing the swab, rotate the swab and repeat the process in the second nostril.

Samples can be taken from both nostrils with the same swab. However, it is not necessary to collect samples from both nostrils if the tip of the swab is saturated with secretion from the first sample collection and a curved nasal septum or blockage causes difficulties in removing the sample from a nostril.

Saliva sample

Avoid consuming food or chewing foreign substances such as food, drink, smoke and chewing gum 30 minutes before collecting the sample.

Place the swab head between the upper and lower molars for 10 seconds to absorb saliva. Press the teeth against the swab head with slight pressure to saturate the swab head with saliva.

Sample preparation

After swab sampling, the swab can be stored in the extraction tube provided with the kit. In addition, the sample can also be stored by immersing the swab head in a vessel containing 2 ml to 3 ml of virus preservation solution, isotonic saline solution, tissue culture solution or phosphate buffer.


Sample transport and storage

Taken samples should be processed as soon as possible and no later than one hour after sampling.

The taken samples must not be stored at 2-8°C for more than 24 hours. At -70°C, the samples can be stored for a longer period of time, but avoid repeated freezing and thawing.

[SAMPLE EXTRACTION]


- Open the closure of the extraction tube and place it in the supplied holder.
- Insert the swab with the nasal swab sample previously taken into the sample tube filled with extraction solution. Roll and rotate the swab at least 5 times while pressing the swab head against the inside and bottom of the sample tube. Move the swab in the sample extraction tube for one minute.
- Remove the sample swab while squeezing the extraction tube with your fingers, squeezing the solution out of the swab head as much as possible. The extracted solution will be used as test sample.
- Close the sample tube with the dripping catch.


[TEST PROCEDURE]

Allow the test and sample to reach room temperature (15 – 30 °C) before performing the test.

- Remove the test kit from the sealed bag.
- Open the cover of the extraction tube and hold it upright while taking 3 drops (approx. 100 µl) of the sample using the dropper pipette and placing it on the well (S) of the test kit. See picture.
- Set the timer to 15 minutes and read the test results after the timer has expired. Do not interpret the test results after 20 minutes.

[INTERPRETATION OF THE RESULTS]

Positive: The colored control line C and test line T appear.

A line appears in the control area (C). A further line appears in the test area (T). A positive result indicated the presence of SARS-CoV-2 nucleocapsid antigen.

Positive results indicate the presence of viral antigens, but a clinical correlation with the patient's medical history and other diagnostic information is required to determine the infection status. Positive results do not rule out a bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Negative: The colored control line C appears, no colored line appears in the test area T.

A line appears in the control area (C). No line appears in the test area (T). Negative results are suspected. Negative test results do not rule out infection and should not be used as the sole basis for treatment or other patient management decisions, including infection control decisions, particularly in the presence of clinical signs and symptoms consistent with COVID-19 or in those who have been in contact with the have been virus. For patient management, it is recommended to confirm these results by a RT-PCR test, if appropriate.

Invalid: The colored control line C does not appear.

Frequent causes for an invalid result are an insufficient sample volume or incorrect process techniques. Check the test procedure and repeat the test with a new kit. If the problem persists, stop using the test kit and contact your dealer or CHASTRU.

[QUALITY CONTROL]

A process control is included in the test. The colored line in the control area (C) is considered to be an internal procedural control. It confirms a sufficient sample volume, a sufficient membrane capillary force and a correct process technique. Control standards are not included with this rapid test kit. It is advisable to carry out positive and negative checks within the framework of good laboratory practice in order to check the test procedure and the correct test performance.

[LIMITATIONS OF THE TEST]

- The COVID-19 rapid antigen test kit is limited to qualitative detection. The intensity of the test line does not necessarily correlate with the concentration of the antigen in the sample. Any coloration in the test line range therefore indicates a positive result.
- Negative results do not exclude SARS-CoV-2 infection and should not be used as the sole criterion for diagnosis.
- A physician must interpret the results in connection with the patient's medical history, physical findings and other diagnostic procedures
- A negative test result is possible if the amount of antigens for the SARS CoV-2 virus in the sample is below the sensitivity threshold value of the test kit or if the target epitope of the virus to which the monoclonal antibody used in the test binds is altered by slight amino acid mutations.
- Proper sampling is required. Incorrect application may lead to inaccurate results. Improper sampling, improper sample storage or repeated freezing and thawing of the samples can lead to inaccurate results.

[FEATURES]

For nasopharyngeal swab samples

Sensitivity, specificity:

The SARS-CoV-2 antigen rapid test was evaluated on nasopharyngeal swab samples from 405 patients and compared with a commercial RT-PCR test. The samples were considered positive if RT-PCR indicated a positive result. The samples were considered negative if RT-PCR indicated a positive result. The RT-PCR comparisons use a nasopharyngeal swab. The result shows the relative sensitivity and specificity.

Method		PCR		Result
CHASTRU Antigen Rapid Test	Result	Positive	Negative	
	Positive	155	1	156
	Negative	6	315	321
Result		161	316	477

Relative sensitivity: 96.34% (95% CI: 92.21%~ 98.65%)

Relative specificity: 99.59% (95% CI: 97.71%~ 99.99%)

For nasal swab sample - anterior nasal

Sensitivity, specificity:

The SARS-CoV-2 Antigen Rapid Test was evaluated on nasopharyngeal swab samples from 405 patients and compared with a commercial RT-PCR test. The samples were considered positive if RT-PCR indicated a positive result. The samples were considered negative if RT-PCR indicated a positive result. The RT-PCR comparisons use a nasopharyngeal swab. The result shows the relative sensitivity and specificity.

Method		PCR		Result
CHASTRU Antigen Rapid Test	Result	Positive	Negative	
	Positive	141	2	143
	Negative	8	412	420
Result		149	414	563

Relative sensitivity: 93.60% (95% CI: 88.20%~ 96.60)

Relative specificity: 99.40% (95% CI: 98.10%~ 99.80%)

For saliva samples

Sensitivity, specificity:

The SARS-CoV-2 Antigen Rapid Test was evaluated on nasopharyngeal swab samples from 405 patients and compared with a commercial RT-PCR test. The samples were considered positive if RT-PCR indicated a positive result. The samples were considered negative if RT-PCR indicated a negative result. The RT-PCR comparisons use a nasopharyngeal swab. The result shows the relative sensitivity and specificity.

Method		PCR		Result
CHASTRU Antigen Rapid Test	Result	Positive	Negative	
	Positive	101	2	103
	Negative	11	423	434
Result		112	425	537

Relative sensitivity: 91.00% (95% CI: 83.80%~ 95.20%)

Relative specificity: 99.40% (95% CI: 98.10%~ 99.80%)

Interfering substances

The following compounds were tested using the SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag), interference could not be detected.

Substance	Concentration	Substance	Concentration
Mucin	2%	Blood	4%
Benzocaine	5mg/mL	Menthol	10mg/ml
Nasal spray (saline)	15%	Phenylephrine	15%
Oxymetazoline	15%	Mupirocin	10mg/ml
Tobramycin	5µg/ml	Zanamivir	5mg/mL
Oseltamivir phosphate	10mg/ml	Ribavirin	5mg/mL
Arbidol	5mg/mL	Dexamethasone	5mg/mL
Fluticasone propionate	5%	Histamine dihydrochloride	10mg/ml
Triamcinolone	10mg/ml		

Limit of detection (analytical sensitivity)

The detection limit of this kit is 50pg/ml recombinant SARS-CoV-2 nucleocapsid protein and 50 TCID50/ml for virus.

Cross-reactivity (analytical specificity)

The cross-reactivity was evaluated by testing 32 commensal and pathogenic microorganisms which may be present in the nasal cavity.

No cross-reactivity with recombinant MERS-CoV-NP protein was observed in tests at a concentration of 50 µg/ml.

No cross-reactivity with the following viruses was observed in tests at a concentration of 1×10^6 PFU/ml: Influenza A (H1N1), Influenza A (H1N1pdm09), Influenza A (H3N2), Influenza B (Yamagata), Influenza B (Victoria), Adenovirus (type 1, 2, 3, 5, 7, 55), Human Metapneumovirus, Parainfluenza virus (type 1, 2, 3, 4), Respiratory Syncytial Virus, Enterovirus, Rhinovirus, Human Coronavirus 229E, Human Coronavirus OC43, Human Coronavirus NL63, Human Coronavirus HKU1.

No cross-reactivity with the following bacteria was observed in tests at a concentration of 1×10^7 PFU/ml: Mycoplasma pneumoniae, Chlamydia pneumoniae, Legionella pneumophila, Haemophilus influenzae, Streptococcus pyogenes (Group A), Streptococcus pneumoniae, Candida albicans, Staphylococcus aureus.


[HOOK EFFECT]

The COVID-19 Antigen Rapid Test kit was tested by a recombinant COVID-19 protein diluted to 100µg/mL and no high-dose hook effect was observed.

Chastru Biotech Limited.

Address: 2/F, Yau Tak Building, 167 Lockhart Road, Wanchai, Hong Kong, China.


Name: SUNGO Europe B.V.

Address: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands

Contact: Yana Zhang (Ms.), Tel/Fax: +31 (0) 2021 11106

Email: yan.zhang@sungoglobal.com

Glossary of symbols

Do not re-use

Temperature limit

Batch code

Use by date

Keep dry

Manufacturer

In vitro diagnostic

Read instructions

Contains sufficient material for <n> tests

Protect from sunlight

Do not use if package is damaged

European Authorized Representation

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