


μCT of lung adenocarcinoma



Monitoring of lung adenocarcinoma progression by ExiTron™ nano 12000

Brett H. Simmons¹, Joseph H. Lee¹, Kush Lalwani¹, Anand Giddabasappa^{1*}, Brittany A. Snider¹, Anthony Wong¹, Patrick B. Lappin¹, Jeetendra Eswaraka¹, Julie L. Kan¹, James G. Christensen¹, and Farbod Shojaei¹

¹ Pfizer Global R&D, 10724 Science Center Dr., La Jolla, San Diego, CA 92121, USA

* Corresponding author (anand.giddabasappa@pfizer.com)

Introduction

The genetically-engineered mouse model, *Kras*^{G12D-LSL} represents a most useful preclinical model system to study tumor growth and progression, constituting a valuable tool to investigate the benefit of new therapies¹. In this study we tested the potency of two drug candidates, PD-0325901 (hereafter PD-901) and PF-04691502 (hereafter PF-502) as single agents or in combination, in suppressing the growth of lung tumors in *Kras*^{G12D-LSL} mice. Previous studies have shown that PF-502 exhibits anti-tumor activity in several preclinical models² and that PD-901 shows potent tumor

growth inhibition in papillary thyroid carcinoma cells³. Also, a combination of both compounds significantly inhibited tumor growth and increased survival compared with either monotherapy in a genetically engineered mouse model of ovarian cancer⁴. Since micro-computed tomography (μCT) has proven to be a reliable method to image lung tumors at different stages of tumorigenesis⁵, the ability of the two drug candidates to suppress lung adenocarcinoma progression was investigated using μCT in combination with the innovative Viscover™ contrast agent ExiTron™ nano 12000.

Materials and methods

Tumor development

Lung tumors were generated in *Kras*^{G12D-LSL} mice (3-4 weeks of age) by intranasally administering an adenovirus vector expressing the Cre recombinase protein (Adeno-Cre)⁶. Standard histological techniques were performed at various timepoints (4, 8, 12, and 16 weeks post-inhalation) to ensure tumor formation and progression.

Therapeutic treatment

Treatment with anti-cancer agents in *Kras*^{G12D-LSL} mice was initiated at 14 weeks after inhalation with Adeno-Cre when lung lesions were at an advanced stage of tumorigenesis (adenocarcinoma). Tumor-bearing mice were treated with anti-cancer compounds, PF-502 (7.5 mg/kg/day), PD-901 (1.5 mg/kg/day) and a combination of PF-502 and PD-901 at the respected dose, for a period of 8 weeks.

μCT Imaging

In vivo μCT of tumor-bearing *Kras*^{G12D-LSL} mice treated with vehicle (blank), PF-502, PD-901, or a combination of PF-502 and PD-901 was performed 4 weeks after treatment initiation at an advanced stage of tumor progression. Imaging was also performed on non-inhaled (control wild type) mice. All mice were anesthetized using isoflurane inhalation and were imaged at 15 min post injection of 100 μL of the CT contrast agent ExiTron™ nano 12000 (Viscover™, nanoPET Pharma GmbH, Berlin, Germany). Upon intravenous injection, ExiTron nano 12000 circulates in the bloodstream




Figure 1: μ CT images obtained 15 min post injection of the CT contrast agent ExiTron nano 12000. **A.** Coronal image of a healthy control wild type mouse. **B.** Transverse images of *Kras*^{G12D-LSL} mice treated with vehicle, PD-901, PF-502, or a combination of PF-502 and PD-901, at 4 weeks after treatment initiation (tumor margins are shown in red).

and accumulates in macrophages, including those in the liver and spleen (reticulo-endothelial system, RES). Due to the enhanced permeation and retention (EPR) effect, the nanoparticulate CT contrast agent can extravasate from fenestrated blood vessels into surrounding tissue of inflamed or tumorous areas and can, thus, assist in the detection of various tumors. All CT images were acquired in standard resolution using respiratory gating (70 kVP, 114 μ A, 300 ms integration time, 41 μ m voxel size).

Results and discussion

To test the efficacy of PF-502 and PD-901 as single agents as well as in combination in the suppression of lung adenocarcinoma, contrast-enhanced *in vivo* μ CT was performed on tumor-bearing mice 4 weeks after initiation of the therapeutic regimens.

Compared with the non-inhaled (control wild type) mice (Fig. 1A), CT images of the vehicle-treated mice clearly showed that the lung tissue is extensively covered by tumor lesions (Fig. 1B) indicating an advanced stage of tumor progression at this timepoint. CT images of mice treated with PD-901, PF-502 and combination therapy exhibited lung tissue with fewer tumors compared to vehicle, indicating that all drug treatments suppressed overall tumor growth. Quantification of the tumor volume clearly showed that while PF-502 partially inhibited tumor growth, PD-901 as single agent or in combination with PF-502 showed superior efficacy (Fig. 2).

These findings were corroborated by histological analysis (data reported elsewhere⁷), implying that the PD-901 compound, as single agent or in combination with PF-502, strongly suppresses lung adenocarcinoma progression. Furthermore, the consistency between data obtained via the two techniques indicates that μ CT using ExiTron nano 12000 is a reliable method to image lung tumors *in vivo*. This is because the contrast agent allows considerable improvement in defining the tumor margins (Fig. 1B, red) due to the reduced contrast intensity of tumor relative to the surrounding blood vessels and heart but increased contrast intensity of tumor relative to air spaces in the lungs. Further work using ExiTron nano 12000 to longitudinally evaluate and quantify lung tumors verifies that this contrast agent allows accurate tumor burden estimation⁸.




Figure 2: Histogram displaying tumor volume (relative to vehicle) in healthy control wild type mice and *Kras*^{G12D-LSL} mice treated with vehicle, PF-502, PD-901, or a combination of PF-502 and PD-901, at 4 weeks after treatment initiation.

Conclusion

In this study, we use contrast-enhanced μ CT to monitor the efficacy of two potential anti-cancer drugs in the suppression of lung adenocarcinoma in $Kras^{G12D-LSL}$ mice. The results show that the applied μ CT method, involving the use of ExiTron nano 12000, enables evaluation of lung tumor progression and treatment strategies in a non-invasive manner, presenting researchers with an efficient tool for anti-cancer drug research and development.

Viscover™ Product	Order No.
ExiTron™ U, 1 x 5 injections	130-095-142
ExiTron™ U, 5 x 5 injections	130-095-143
ExiTron™ V, 1 x 5 injections	130-095-283
ExiTron™ V, 5 x 5 injections	130-095-284
ExiTron™ P, 1 x 5 injections	130-095-144
ExiTron™ P, 5 x 5 injections	130-095-145
ExiTron™ nano 6000, 1 x 5 injections	130-095-146
ExiTron™ nano 6000, 5 x 5 injections	130-095-147
ExiTron™ nano 12000, 1 x 5 injections	130-095-698
ExiTron™ nano 12000, 5 x 5 injections	130-095-700
ExiTron™ MyoC 8000, 1 x 5 injections	130-095-701
ExiTron™ MyoC 8000, 5 x 5 injections	130-095-702

References

1. Jackson, E.L. *et al.* (2001) Analysis of lung tumor initiation and progression using conditional expression of oncogenic *K-ras*. *Genes Dev.* 15(24): 3243–3248.
2. Yuan, J. *et al.* (2011) PF-04691502, a potent and selective oral inhibitor of PI3K and mTOR kinases with antitumor activity. *Mol Cancer Ther.* 10(11): 2189-2199.
3. Henderson, Y.C. *et al.* (2010) MEK inhibitor PD0325901 significantly reduces the growth of papillary thyroid carcinoma cells *in vitro* and *in vivo*. *Mol Cancer Ther.* 9(7): 1968–1976.
4. Kinross, K.M. *et al.* (2011) *In Vivo* Activity of Combined PI3K/mTOR and MEK Inhibition in a *Kras*^{G12D};*Pten* Deletion Mouse Model of Ovarian Cancer. *Mol Cancer Ther.* 10(8): 1440–1449.
5. Ritman, E.L. (2011) Current status of developments and applications of micro-CT. *Annu Rev Biomed Eng.* 13: 531–552.
6. DuPage, M. *et al.* (2009) Conditional mouse lung cancer models using adenoviral or lentiviral delivery of Cre recombinase. *Nat Protoc.* 4(7): 1064–1072.
7. Simmons, B.H. *et al.* (2012) Combination of a MEK inhibitor at sub-MTD with a PI3K/mTOR inhibitor significantly suppresses growth of lung adenocarcinoma tumors in $Kras^{G12D-LSL}$ mice. *Cancer Chemother Pharmacol.* 70(2): 213-220.
8. Lalwani, K. *et al.* (2013) Contrast agents for quantitative microCT of lung tumors in mice. *Comp Med.* 63(6): 482-490.

VISCOVER

powered by

nanOPET

nanoPET Pharma GmbH

Robert-Koch-Platz 4

10115 Berlin, Germany

Phone: +49 30 890 49 74 - 0

Fax: +49 30 890 49 74 - 99

Email: imaging@viscover.berlin

nanoPET Pharma GmbH provides products and services worldwide. Visit **www.viscover.berlin** to find your nearest sales contact.

Unless otherwise specifically indicated, *nanoPET Pharma GmbH* products and services are for research use only and not for therapeutic or diagnostic use. ExiTron and Viscover are trademarks of *nanoPET Pharma GmbH*. All other trademarks mentioned in this document are the property of their respective owners and are used for identification purposes only.

Copyright © 2018 *nanoPET Pharma GmbH*. All rights reserved.