





Uroflex® Urinalkondom Schablone zum Ausschneiden



Mit dieser Schablone finden Sie ganz einfach die richtige Größe für Ihr Urinalkondom.

Uroflex® Urinalkondome Typ I bis III

REF	Typ	Gesamtlänge mm	Schaftlänge mm	Klebelänge mm	Innen Ø mm	PZN	VE in Stk.
3A0125	 Typ I	145	80	45	25	17414651	30
3A0129		145	80	45	29	17414668	30
3A0132		145	80	45	32	17414674	30
3A0136		145	80	45	36	17414680	30
3A0141		145	80	45	41	17414697	30
3A0225	 Typ II	115	50	40	25	17414705	30
3A0229		115	50	40	29	17414711	30
3A0232		115	50	40	32	17414728	30
3A0236		115	50	40	36	17414734	30
3A0241		115	50	40	41	17414740	30
3A0325	 Typ III	145	80	70	25	17414757	30
3A0329		145	80	70	29	17414763	30
3A0332		145	80	70	32	17414786	30
3A0336		145	80	70	36	17414792	30
3A0341		145	80	70	41	17414800	30
(± 5mm)							